

# Recombinant Human DNAJC12

Catalog No: #GP13278

Package Size: #GP13278-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human DNAJC12
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-198 amino acids of human DNAJC12
Target Name	DnaJ heat shock protein family (Hsp40) member C12
Other Names	JDP1; HPANBH4
Accession No.	Swissprot:Q9UKB3 Gene Accession:BC017018
Uniprot	Q9UKB3
GeneID	56521;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only