

Recombinant Human RTN1

Catalog No: #GP13280

Package Size: #GP13280-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RTN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 589-776 amino acids of human RTN1
Target Name	reticulon 1
Other Names	NSP
Accession No.	Swissprot:Q16799 Gene Accession:BC000314
Uniprot	Q16799
GeneID	6252;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the family of reticulon encoding genes. Reticulons are associated with the endoplasmic reticulum, and are involved in neuroendocrine secretion or in membrane trafficking in neuroendocrine cells. This gene is considered to be a specific marker for neurological diseases and cancer, and is a potential molecular target for therapy. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only