

Recombinant Human ARL4D

Catalog No: #GP13281



Package Size: #GP13281-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARL4D
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-201 amino acids of human ARL4D
Target Name	ADP ribosylation factor like GTPase 4D
Other Names	ARL6; ARF4L
Accession No.	Swissprot:P49703Gene Accession:BC000043
Uniprot	P49703
GeneID	379;
Storage	-20~-80°C, pH 7.6 PBS

Background

ADP-ribosylation factor 4D is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL4D is closely similar to ARL4A and ARL4C and each has a nuclear localization signal and an unusually high guanine nucleotide exchange rate. This protein may play a role in membrane-associated intracellular trafficking. Mutations in this gene have been associated with Bardet-Biedl syndrome (BBS).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only