

Recombinant Human RAB1A

Catalog No: #GP13284



Package Size: #GP13284-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human RAB1A |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 6-205 amino acids of human RAB1A |
| Target Name | RAB1A, member RAS oncogene family |
| Other Names | RAB1; YPT1 |
| Accession No. | Swissprot:P62820Gene Accession:BC000905 |
| Uniprot | P62820 |
| GeneID | 5861; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a member of the Ras superfamily of GTPases. Members of the gene family cycle between inactive GDP-bound and active GTP-bound forms. This small GTPase controls vesicle traffic from the endoplasmic reticulum to the Golgi apparatus. Multiple alternatively spliced transcript variants have been identified for this gene which encode different protein isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only