

Recombinant Human HS1BP3

Catalog No: #GP13286



Package Size: #GP13286-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HS1BP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 193-392 amino acids of human HS1BP3
Target Name	HCLS1 binding protein 3
Other Names	ETM2; HS1-BP3
Accession No.	Swissprot:Q53T59Gene Accession:BC050636
Uniprot	Q53T59
GeneID	64342;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only