

Recombinant Human PAFAH1B2

Catalog No: #GP13298

Package Size: #GP13298-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAFAH1B2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-229 amino acids of human PAFAH1B2
Target Name	platelet activating factor acetylhydrolase 1b catalytic subunit 2
Other Names	HEL-S-303
Accession No.	Swissprot:P68402Gene Accession:BC000398
Uniprot	P68402
GeneID	5049;
Storage	-20~-80°C, pH 7.6 PBS

Background

Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only