

Recombinant Human THTPA

Catalog No: #GP13299

Package Size: #GP13299-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human THTPA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-230 amino acids of human THTPA
Target Name	thiamine triphosphatase
Other Names	THTP; THTPASE
Accession No.	Swissprot:Q9BU02Gene Accession:BC002984
Uniprot	Q9BU02
GeneID	79178;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an enzyme which catalyzes the biosynthesis of thiamine diphosphate (vitamin B1) by hydrolysis of thiamine triphosphate. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only