

Recombinant Human PDCL2

Catalog No: #GP13302



Package Size: #GP13302-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDCL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-239 amino acids of human PDCL2
Target Name	phosducin like 2
Other Names	GCPHLP
Accession No.	Swissprot:Q8N4E4Gene Accession:BC034431
Uniprot	Q8N4E4
GeneID	132954;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the phosducin-like protein family and is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin. Members of the phosducin-like protein family have been shown to bind to the beta-gamma subunits of G proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only