

Recombinant Human SULT4A1

Catalog No: #GP13303



Package Size: #GP13303-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SULT4A1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-260 amino acids of human SULT4A1
Target Name	sulfotransferase family 4A member 1
Other Names	NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1
Accession No.	Swissprot:Q9BR01 Gene Accession:BC028171
Uniprot	Q9BR01
GeneID	25830;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only