

Recombinant Human HCCS

Catalog No: #GP13305



Package Size: #GP13305-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HCCS
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-268 amino acids of human HCCS
Target Name	holocytochrome c synthase
Other Names	MLS; CCHL; MCOPS7; LSDMCA1
Accession No.	Swissprot:P53701Gene Accession:BC001691
Uniprot	P53701
GeneID	3052;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is an enzyme that covalently links a heme group to the apoprotein of cytochrome c. Defects in this gene are a cause of microphthalmia syndromic type 7 (MCOPS7). Three transcript variants encoding the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only