

Recombinant Human RETREG1

Catalog No: #GP13320



Package Size: #GP13320-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RETREG1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 357-497 amino acids of human RETREG1
Target Name	reticulophagy regulator 1
Other Names	JK1; JK-1; FAM134B
Accession No.	Swissprot:Q9H6L5Gene Accession:BC030517
Uniprot	Q9H6L5
GeneID	54463;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a cis-Golgi transmembrane protein that may be necessary for the long-term survival of nociceptive and autonomic ganglion neurons. Mutations in this gene are a cause of hereditary sensory and autonomic neuropathy type IIB (HSAN IIB), and this gene may also play a role in susceptibility to vascular dementia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only