

Recombinant Human ASTE1

Catalog No: #GP13323



Package Size: #GP13323-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ASTE1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 480-679 amino acids of human ASTE1
Target Name	asteroid homolog 1
Other Names	HT001
Accession No.	Swissprot:Q2TB18Gene Accession:BC030274
Uniprot	Q2TB18
GeneID	28990;
Storage	-20~-80°C, pH 7.6 PBS

Background

Possible role in EGF receptor signaling.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only