

Recombinant Human TEAD3

Catalog No: #GP13325



Package Size: #GP13325-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TEAD3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 236-435 amino acids of human TEAD3
Target Name	TEA domain transcription factor 3
Other Names	TEF5; TEAD5; TEF-5; DTEF-1; ETFR-1; TEAD-3
Accession No.	Swissprot:Q99594Gene Accession:BC027877
Uniprot	Q99594
GeneID	7005;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only