

Recombinant Human CYR1

Catalog No: #GP13329



Package Size: #GP13329-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CYR1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 30-154 amino acids of human CYR1
Target Name	cysteine and tyrosine rich 1
Other Names	C21orf95
Accession No.	Swissprot:Q96J86Gene Accession:BC036761
Uniprot	Q96J86
GeneID	116159;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein is unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only