

Recombinant Human KCTD13

Catalog No: #GP13330



Package Size: #GP13330-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KCTD13
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 130-329 amino acids of human KCTD13
Target Name	potassium channel tetramerization domain containing 13
Other Names	PDIP1; FKSG86; BACURD1; POLDIP1; hBACURD1
Accession No.	Swissprot:Q8WZ19Gene Accession:BC036228
Uniprot	Q8WZ19
GeneID	253980;
Storage	-20~-80°C, pH 7.6 PBS

Background

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex involved in regulation of cytoskeleton structure. The BCR(BACURD1) E3 ubiquitin ligase complex mediates the ubiquitination of RHOA, leading to its degradation by the proteasome, thereby regulating the actin cytoskeleton and cell migration.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only