

Recombinant Human PDS5A

Catalog No: #GP13331



Package Size: #GP13331-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDS5A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1138-1337 amino acids of human PDS5A
Target Name	PDS5 cohesin associated factor A
Other Names	PIG54; SCC112; SCC-112
Accession No.	Swissprot:Q29RF7 Gene Accession:BC009650
Uniprot	Q29RF7
GeneID	23244;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene binds to the cohesin complex and associates with chromatin through most of the cell cycle. The encoded protein may play a role in regulating sister chromatid cohesion during mitosis. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only