

Recombinant Human C19orf57

Catalog No: #GP13332



Package Size: #GP13332-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human C19orf57 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 471-658 amino acids of human C19orf57 |
| Target Name | chromosome 19 open reading frame 57 |
| Accession No. | Swissprot:Q0VDD7Gene Accession:BC012945 |
| Uniprot | Q0VDD7 |
| GeneID | 79173; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only