

# Recombinant Human IQCD

Catalog No: #GP13341



Package Size: #GP13341-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human IQCD
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human IQCD
Target Name	IQ motif containing D
Other Names	DRC10; CFAP84; 4933433C09Rik
Accession No.	Swissprot:Q96DY2Gene Accession:BC013151
Uniprot	Q96DY2
GeneID	115811;
Storage	-20~-80°C, pH 7.6 PBS

## Background

IQCD contains 1 IQ domain. The function of the IQCD protein is not known.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only