

Recombinant Human LRG1

Catalog No: #GP13342



Package Size: #GP13342-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LRG1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 148-347 amino acids of human LRG1
Target Name	leucine rich alpha-2-glycoprotein 1
Other Names	LRG; HMFT1766
Accession No.	Swissprot:P02750Gene Accession:BC034389
Uniprot	P02750
GeneID	116844;
Storage	-20~-80°C, pH 7.6 PBS

Background

The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002 [PubMed 12223515]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only