

Recombinant Human TRUB1

Catalog No: #GP13343



Package Size: #GP13343-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TRUB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 150-349 amino acids of human TRUB1
Target Name	TruB pseudouridine synthase family member 1
Other Names	PUS4
Accession No.	Swissprot:Q8WWH5Gene Accession:BC030601
Uniprot	Q8WWH5
GeneID	142940;
Storage	-20~-80°C, pH 7.6 PBS

Background

Pseudouridine is an abundant component of rRNAs and tRNAs and is enzymatically generated by isomerization of uridine by pseudouridine synthase (Zucchini et al., 2003 [PubMed 12736709]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only