

## Recombinant Human MKX

Catalog No: #GP13345



Package Size: #GP13345-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MKX
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 196-352 amino acids of human MKX
Target Name	mohawk homeobox
Other Names	IFRX; IRXL1; C10orf48
Accession No.	Swissprot:Q8IYA7Gene Accession:BC036207
Uniprot	Q8IYA7
GeneID	283078;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is an IRX family-related homeobox protein that may play a role in cell adhesion. Studies in mice suggest that this protein may be a regulator of tendon development. Two transcript variants encoding the same protein have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only