

# Recombinant Human ELAC1

Catalog No: #GP13351



Package Size: #GP13351-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ELAC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 164-363 amino acids of human ELAC1
Target Name	elaC ribonuclease Z 1
Other Names	D29
Accession No.	Swissprot:Q9H777Gene Accession:BC014624
Uniprot	Q9H777
GeneID	55520;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Zinc phosphodiesterase, which displays some tRNA 3'-processing endonuclease activity. Probably involved in tRNA maturation, by removing a 3'-trailer from precursor tRNA.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only