

Recombinant Human RABEPK

Catalog No: #GP13357



Package Size: #GP13357-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RABEPK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 173-372 amino acids of human RABEPK
Target Name	Rab9 effector protein with kelch motifs
Other Names	p40; RAB9P40; bA65N13.1
Accession No.	Swissprot:Q7Z6M1 Gene Accession:BC000503
Uniprot	Q7Z6M1
GeneID	10244;
Storage	-20~-80°C, pH 7.6 PBS

Background

Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only