

Recombinant Human SAMSN1

Catalog No: #GP13358



Package Size: #GP13358-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SAMSN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 174-373 amino acids of human SAMSN1
Target Name	SAM domain, SH3 domain and nuclear localization signals 1
Other Names	SLy2; HACS1; NASH1; SASH2; SH3D6B
Accession No.	Swissprot:Q9NSI8Gene Accession:BC029112
Uniprot	Q9NSI8
GeneID	64092;
Storage	-20~-80°C, pH 7.6 PBS

Background

SAMSN1 is a member of a novel gene family of putative adaptors and scaffold proteins containing SH3 and SAM (sterile alpha motif) domains (Claudio et al., 2001 [PubMed 11536050]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only