

# Recombinant Human MAGEC2

Catalog No: #GP13360



Package Size: #GP13360-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MAGEC2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 174-373 amino acids of human MAGEC2
Target Name	MAGE family member C2
Other Names	CT10; HCA587; MAGEE1
Accession No.	Swissprot:Q9UBF1 Gene Accession:BC013318
Uniprot	Q9UBF1
GeneID	51438;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the MAGEC gene family. It is not expressed in normal tissues, except for testis, and is expressed in tumors of various histological types. This gene and the other MAGEC genes are clustered on chromosome Xq26-q27.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only