

Recombinant Human LDB1

Catalog No: #GP13361



Package Size: #GP13361-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LDB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 212-411 amino acids of human LDB1
Target Name	LIM domain binding 1
Other Names	NLI; CLIM2; LDB-1; CLIM-2
Accession No.	Swissprot:Q86U70Gene Accession:BC000482
Uniprot	Q86U70
GeneID	8861;
Storage	-20~-80°C, pH 7.6 PBS

Background

Binds to the LIM domain of a wide variety of LIM domain-containing transcription factors. May regulate the transcriptional activity of LIM-containing proteins by determining specific partner interactions. Plays a role in the development of interneurons and motor neurons in cooperation with LHX3 and ISL1. Acts synergistically with LHX1/LIM1 in axis formation and activation of gene expression. Acts with LMO2 in the regulation of red blood cell development, maintaining erythroid precursors in an immature state (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only