

Recombinant Human LDHAL6B

Catalog No: #GP13364



Package Size: #GP13364-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LDHAL6B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 182-381 amino acids of human LDHAL6B
Target Name	lactate dehydrogenase A like 6B
Other Names	LDHL; LDH6B; LDHAL6
Accession No.	Swissprot:Q9BYZ2Gene Accession:BC022034
Uniprot	Q9BYZ2
GeneID	92483;
Storage	-20~-80°C, pH 7.6 PBS

Background

LDHAL6B belongs to the LDH/MDH superfamily. LDH family. The specific function of the protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only