

## Recombinant Human PCYT2

Catalog No: #GP13369



Package Size: #GP13369-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PCYT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 190-389 amino acids of human PCYT2
Target Name	phosphate cytidyltransferase 2, ethanolamine
Other Names	ET
Accession No.	Swissprot:Q99447Gene Accession:BC000351
Uniprot	Q99447
GeneID	5833;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an enzyme that catalyzes the formation of CDP-ethanolamine from CTP and phosphoethanolamine in the Kennedy pathway of phospholipid synthesis. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only