

Recombinant Human TENT5C

Catalog No: #GP13371



Package Size: #GP13371-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TENT5C
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 192-391 amino acids of human TENT5C
Target Name	terminal nucleotidyltransferase 5C
Other Names	FAM46C
Accession No.	Swissprot:Q5VWP2Gene Accession:BC036516
Uniprot	Q5VWP2
GeneID	54855;
Storage	-20~-80°C, pH 7.6 PBS

Background

Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. (Microbial infection) Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only