

Recombinant Human TUBD1

Catalog No: #GP13374



Package Size: #GP13374-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TUBD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human TUBD1
Target Name	tubulin delta 1
Other Names	TUBD
Accession No.	Swissprot:Q9UJT1 Gene Accession:BC000258
Uniprot	Q9UJT1
GeneID	51174;
Storage	-20~-80°C, pH 7.6 PBS

Background

In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only