

Recombinant Human UBAC1

Catalog No: #GP13377



Package Size: #GP13377-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBAC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 206-405 amino acids of human UBAC1
Target Name	UBA domain containing 1
Other Names	KPC2; GBDR1; UBADC1
Accession No.	Swissprot:Q9BSL1 Gene Accession:BC004967
Uniprot	Q9BSL1
GeneID	10422;
Storage	-20~-80°C, pH 7.6 PBS

Background

Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only