

Recombinant Human UXS1

Catalog No: #GP13382



Package Size: #GP13382-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UXS1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 195-394 amino acids of human UXS1
Target Name	UDP-glucuronate decarboxylase 1
Other Names	UGD; SDR6E1
Accession No.	Swissprot:Q8NBZ7 Gene Accession:BC009819
Uniprot	Q8NBZ7
GeneID	80146;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an enzyme found in the perinuclear Golgi which catalyzes the synthesis of UDP-xylose used in glycosaminoglycan (GAG) synthesis on proteoglycans. The GAG chains are covalently attached to proteoglycans which participate in signaling pathways during development. Multiple transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only