

Recombinant Human AK9

Catalog No: #GP13384



Package Size: #GP13384-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AK9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 222-421 amino acids of human AK9
Target Name	adenylate kinase 9
Other Names	AK 9; AKD1; AKD2; C6orf199; C6orf224; dJ70A9.1
Accession No.	Swissprot:Q5TCS8Gene Accession:BC022031
Uniprot	Q5TCS8
GeneID	221264;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene catalyzes the interconversion of nucleosides, possessing both nucleoside monophosphate and diphosphate kinase activities. The encoded protein uses these interconversions to maintain nucleoside homeostasis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only