

Recombinant Human ECSIT

Catalog No: #GP13392



Package Size: #GP13392-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ECSIT
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 232-431 amino acids of human ECSIT
Target Name	ECSIT signalling integrator
Other Names	SITPEC
Accession No.	Swissprot:Q9BQ95Gene Accession:BC000193
Uniprot	Q9BQ95
GeneID	51295;
Storage	-20~-80°C, pH 7.6 PBS

Background

Adapter protein of the Toll-like and IL-1 receptor signaling pathway that is involved in the activation of NF-kappa-B via MAP3K1. Promotes proteolytic activation of MAP3K1. Involved in the BMP signaling pathway. Required for normal embryonic development (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only