

# Recombinant Human RCOR3

Catalog No: #GP13394



Package Size: #GP13394-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RCOR3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 59-258 amino acids of human RCOR3
Target Name	REST corepressor 3
Accession No.	Swissprot:Q9P2K3Gene Accession:BC031608
Uniprot	Q9P2K3
GeneID	55758;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May act as a component of a corepressor complex that represses transcription.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only