

Recombinant Human SIGLEC6

Catalog No: #GP13395



Package Size: #GP13395-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SIGLEC6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 254-453 amino acids of human SIGLEC6
Target Name	sialic acid binding Ig like lectin 6
Other Names	CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327
Accession No.	Swissprot:O43699Gene Accession:BC035359
Uniprot	O43699
GeneID	946;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the SIGLEC (sialic acid binding immunoglobulin-like lectin) family of proteins. The encoded transmembrane receptor binds sialyl-TN glycans and leptin. Placental expression of the encoded protein is upregulated in preeclampsia.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only