

Recombinant Human ASTN2

Catalog No: #GP13396



Package Size: #GP13396-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ASTN2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1140-1339 amino acids of human ASTN2
Target Name	astrotactin 2
Other Names	bA67K19.1
Accession No.	Swissprot:O75129Gene Accession:BC029272
Uniprot	O75129
GeneID	23245;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is expressed in the brain and may function in neuronal migration, based on functional studies of the related astrotactin 1 gene in human and mouse. A deletion at this locus has been associated with schizophrenia. Multiple transcript variants encoding different proteins have been found for this locus.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only