

Recombinant Human TUBB2A

Catalog No: #GP13398



Package Size: #GP13398-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TUBB2A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 246-445 amino acids of human TUBB2A
Target Name	tubulin beta 2A class IIa
Other Names	TUBB; TUBB2; CDCBM5
Accession No.	Swissprot:Q13885Gene Accession:BC001194
Uniprot	Q13885
GeneID	7280;
Storage	-20~-80°C, pH 7.6 PBS

Background

Microtubules, key participants in processes such as mitosis and intracellular transport, are composed of heterodimers of alpha- and beta-tubulins. The protein encoded by this gene is a beta-tubulin. Defects in this gene are associated with complex cortical dysplasia with other brain malformations-5. Two transcript variants encoding distinct isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only