

Recombinant Human AKAP17A

Catalog No: #GP13399



Package Size: #GP13399-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AKAP17A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human AKAP17A
Target Name	A-kinase anchoring protein 17A
Other Names	XE7; 721P; XE7Y; CCDC133; CXYorf3; PRKA17A; SFRS17A; AKAP-17A; DXYS155E
Accession No.	Swissprot:Q02040Gene Accession:BC028151
Uniprot	Q02040
GeneID	8227;
Storage	-20~-80°C, pH 7.6 PBS

Background

This locus encodes a protein kinase A anchoring protein. The encoded protein is part of the spliceosome complex and is involved in the regulation of alternate splicing in some mRNA precursors. Alternatively spliced transcript variants have been identified for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only