

Recombinant Human NADK

Catalog No: #GP13400



Package Size: #GP13400-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NADK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 247-446 amino acids of human NADK
Target Name	NAD kinase
Other Names	dJ283E3.1
Accession No.	Swissprot:O95544Gene Accession:BC001709
Uniprot	O95544
GeneID	65220;
Storage	-20~-80°C, pH 7.6 PBS

Background

NADK catalyzes the transfer of a phosphate group from ATP to NAD to generate NADP, which in its reduced form acts as an electron donor for biosynthetic reactions (Lerner et al., 2001 [PubMed 11594753]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only