

## Recombinant Human INPP5K

Catalog No: #GP13401



Package Size: #GP13401-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human INPP5K
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 249-448 amino acids of human INPP5K
Target Name	inositol polyphosphate-5-phosphatase K
Other Names	PPS; SKIP; MDCCAID
Accession No.	Swissprot:Q9BT40Gene Accession:BC004362
Uniprot	Q9BT40
GeneID	51763;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only