

# Recombinant Human GALK2

Catalog No: #GP13406



Package Size: #GP13406-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GALK2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 259-458 amino acids of human GALK2
Target Name	galactokinase 2
Other Names	GK2
Accession No.	Swissprot:Q01415Gene Accession:BC005141
Uniprot	Q01415
GeneID	2585;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. The encoded protein is a member of the GHMP kinase family. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only