

# Recombinant Human ZNF256

Catalog No: #GP13414



Package Size: #GP13414-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZNF256
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 428-627 amino acids of human ZNF256
Target Name	zinc finger protein 256
Other Names	BMZF3; BMZF-3
Accession No.	Swissprot:Q9Y2P7Gene Accession:BC001438
Uniprot	Q9Y2P7
GeneID	10172;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Transcriptional repressor that plays a role in cell proliferation. Requires TRIM28 for its activity.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only