

Recombinant Human GPKOW

Catalog No: #GP13415



Package Size: #GP13415-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GPKOW
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 277-476 amino acids of human GPKOW
Target Name	G-patch domain and KOW motifs
Other Names	T54; Spp2; GPATC5; GPATCH5
Accession No.	Swissprot:Q92917Gene Accession:BC000397
Uniprot	Q92917
GeneID	27238;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a putative RNA-binding protein containing G-patch and KOW (Kyprides, Ouzounis, Woese) domains. The encoded protein interacts directly with protein kinase A and protein kinase X and is also found associated with the spliceosome.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only