

## Recombinant Human FDXR

Catalog No: #GP13418



Package Size: #GP13418-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FDXR
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 292-491 amino acids of human FDXR
Target Name	ferredoxin reductase
Other Names	ADXR; ANOA
Accession No.	Swissprot:P22570Gene Accession:BC002960
Uniprot	P22570
GeneID	2232;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a mitochondrial flavoprotein that initiates electron transport for cytochromes P450 receiving electrons from NADPH. Multiple alternatively spliced transcript variants have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only