

# Recombinant Human ARIH2

Catalog No: #GP13421



Package Size: #GP13421-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARIH2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 294-493 amino acids of human ARIH2
Target Name	ariadne RBR E3 ubiquitin protein ligase 2
Other Names	ARI2; TRIAD1
Accession No.	Swissprot:O95376Gene Accession:BC000422
Uniprot	O95376
GeneID	10425;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is an E3 ubiquitin-protein ligase that polyubiquitinates some proteins, tagging them for degradation. The encoded protein upregulates p53 in some cancer cells and may inhibit myelopoiesis. Several transcript variants encoding different isoforms have been found for this gene, although the full-length nature of some of them have not been determined yet.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only