

Recombinant Human AHCYL1

Catalog No: #GP13424



Package Size: #GP13424-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AHCYL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 331-530 amino acids of human AHCYL1
Target Name	adenosylhomocysteinase like 1
Other Names	DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA
Accession No.	Swissprot:O43865Gene Accession:BC007576
Uniprot	O43865
GeneID	10768;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene interacts with inositol 1,4,5-trisphosphate receptor, type 1 and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding two different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only