

# Recombinant Human ZNF559

Catalog No: #GP13428



Package Size: #GP13428-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZNF559
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 339-538 amino acids of human ZNF559
Target Name	zinc finger protein 559
Other Names	NBLA00121
Accession No.	Swissprot:Q9BR84Gene Accession:BC006436
Uniprot	Q9BR84
GeneID	84527;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May be involved in transcriptional regulation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only