

Recombinant Human RIOK1

Catalog No: #GP13431



Package Size: #GP13431-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RIOK1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 369-568 amino acids of human RIOK1
Target Name	RIO kinase 1
Other Names	RIO1; AD034; RRP10; bA288G3.1
Accession No.	Swissprot:Q9BRS2Gene Accession:BC006104
Uniprot	Q9BRS2
GeneID	83732;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene competes with pICln for inclusion in the protein arginine methyltransferase 5 complex. This complex targets substrates for dimethylation. The encoded protein is essential for the last steps in the maturation of 40S subunits.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only