

Recombinant Human PAK1IP1

Catalog No: #GP13433



Package Size: #GP13433-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAK1IP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 409-608 amino acids of human DENND6A
Target Name	DENN domain containing 6A
Other Names	AF11A; FAM116A
Accession No.	Swissprot:Q8IWF6Gene Accession:BC040291
Uniprot	Q8IWF6
GeneID	201627;
Storage	-20~-80°C, pH 7.6 PBS

Background

Guanine nucleotide exchange factor (GEF) for RAB14. Component of an endocytic recycling pathway that is required for the control of ADAM10 transport, shedding of N-cadherin/CDH2 by ADAM9 or ADAM10 and regulation of cell-cell junctions. Required for RAB14 recruitment to recycling endosomes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only